



**BUBBLE CPAP & NON-INVASIVE  
RESPIRATORY MANAGEMENT OF THE NEWBORN**  
*Conference & Workshop*

# Bubble CPAP from the Respiratory Therapy Perspective

Christopher Thornton, RRT  
Respiratory Therapist  
Anne Arundel Medical Center  
Annapolis, MD

# **BUBBLE CPAP & NON-INVASIVE RESPIRATORY MANAGEMENT OF THE NEWBORN**

*Conference & Workshop*

---



## **Disclosures**

Off-Label Usage: None

Interests: None

# Building Bubble CPAP

## Objectives:

1. Demonstrate Building of Bubble CPAP
2. Guide thru 'Hands-on' building of CPAP
3. Explain Gas flow thru Bubble CPAP
4. Discuss Troubleshooting
5. Parts list for Bubble CPAP
6. Questions

# Building Bubble CPAP



# Building Bubble CPAP

## 5 Building Steps

1. Blender
2. Humidifier
3. Inspiratory & Expiratory tubing
4. Bubble bottle
5. Test

# Building Bubble CPAP

- 1) Blender: Controls
  - a.  $FIO_2$
  - b. Flow rate



# Building Bubble CPAP

2. Humidifier assembly
  - a. Connects to Blender
  - b. Access for O<sub>2</sub> Analyzer
  - c. Filters gas source
  - d. Humidifies Gas



# Building Bubble CPAP

## 2. Humidifier Assembly




5 Pieces Required



# Building Bubble CPAP


## 2. Humidifier Assembly



1) Oxygen Barb  $\frac{3}{4}$ " O.D.

# Building Bubble CPAP


## 2. Humidifier Assembly



2) Adapter, Trach T 22mm

# Building Bubble CPAP


## 2. Humidifier Assembly



## 3) Bacterial/ Viral Filter

# Building Bubble CPAP

## 2. Humidifier Assembly



4) Airlife Cuff Connector  
22 mm I.D.

# Building Bubble CPAP

## 2. Humidifier Assembly



5) Humidifier Chamber

# Building Bubble CPAP

## 2. Humidifier Assembly

- a. Attach Oxygen Barb to Adapter, Trach T
- b. Attach Filter to above
- c. Attach Airlife Cuff Connector to above
- d. Attach above to Humidifier Chamber

Find pieces & Assemble



# Building Bubble CPAP

## 3. Inspiratory & Expiratory Tubing

- a. Blue tubing moves Inspiratory Gas to baby
- b. White tubing moves Expiratory Gas to Bubble Bottle
- c. Blue tubing by convention is Inspiratory
- d. White tubing by convention is Expiratory

F&P 750 Circuit



F&P 850 Circuit



# Building Bubble CPAP

## 4. Bubble Bottle

- a. Connects to Expiratory Tubing
- b. Holds Acetic Acid solution
- c. Capillary tube depth provides CPAP





# Building Bubble CPAP


## 4. Bubble Bottle



5 Pieces Required

# Building Bubble CPAP

## 4. Bubble Bottle



1) Measuring tape (cut from 0 to 8cm)

# Building Bubble CPAP


## 4. Bubble Bottle



## 2) Capillary Tube

# Building Bubble CPAP


## 4. Bubble Bottle



3) Mask Intubation Adapter(15mm X 22mm)

# Building Bubble CPAP

## 4. Bubble Bottle



4) Silicone Flex-adapter

# Building Bubble CPAP

## 4. Bubble Bottle



## 5) 0.25% Acetic Acid Bottle

# Building Bubble CPAP

## 4. Bubble Bottle



- a. Attach Measuring Tape 8cm at bottom with Zero up
- b. Attach Capillary tube to Mask Intubation adapter
- c. Thread bottom of Capillary tube thru Silicone Flex-adapter and push  $\frac{1}{2}$  way up and over Mask Intubation Adapter
- d. Place Capillary Assembly in Bottle top
- e. Empty Acetic Acid to Zero
- f. Cut Gas escape hole in Bottle

Find Pieces and Assemble

# Building Bubble CPAP

## Final build Bubble CPAP System



1) Blender



# Building Bubble CPAP

Final build Bubble CPAP System



1) Blender



Attach Blender  
to Humidifier Assembly  
with Supply Tubing



2) Humidifier  
Assembly

# Building Bubble CPAP

## Final build Bubble CPAP System



1) Blender



2) Humidifier Assembly

Supply Tubing



**BABY**

3) Inspiratory Tubing



Attach  
Inspiratory  
Tubing(Blue)  
To  
Humidifier  
Assembly

# Building Bubble CPAP

Final build Bubble CPAP System



1) Blender



2) Humidifier Assembly

Supply Tubing

**BABY**

3) Inspiratory Tubing



3) Expiratory Tubing

4) Bubble Bottle Assembly



Attach  
Expiratory  
Tubing(White)  
To  
Bubble Bottle  
Assembly

# Building Bubble CPAP

## Gas Flow thru the Bubble CPAP System



1) Blender



2) Humidifier Assembly

**BABY**

**Atmosphere**


Supply Tubing

3) Inspiratory Tubing



3) Expiratory Tubing

4) Bubble Bottle Assembly



# Building Bubble CPAP

**To Test:** Turn on Flow Meter if it Bubbles.....



**Its Working!!**

# Building Bubble CPAP

Troubleshooting:

“Not Bubbling” equals a LEAK

- a. #1 leak is the babies mouth
- b. Connections are not tight
- c. CPAP Prongs are too small
- d. Check tubing connections starting at Blender to Bubble Bottle for gas escaping

# MOC Question #1

Gas flow through the Bubble CPAP system could best be described by which of the following?

- A) Gas source to heater to Bubble CPAP device to inspiratory limb to baby to expiratory limb to to atmosphere
- B) Bubble CPAP device to heater to inspiratory limb to baby to expiratory limb to atmosphere
- C) Gas source to heater to inspiratory limb to baby to expiratory limb to Bubble CPAP device to atmosphere
- D) None of the above

# MOC Answer #1

Gas flow through the Bubble CPAP system could best be described by which of the following?

- A) Gas source to heater to Bubble CPAP device to inspiratory limb to baby to expiratory limb to to atmosphere
- B) Bubble CPAP device to heater to inspiratory limb to baby to expiratory limb to atmosphere
- C) **Gas source to heater to inspiratory limb to baby to expiratory limb to Bubble CPAP device to atmosphere**
- D) None of the above



# Building Bubble CPAP

## Gas Flow thru the Bubble CPAP System



1) Blender



2) Humidifier Assembly

**BABY**

**Atmosphere**

Supply Tubing

3) Inspiratory Tubing



3) Expiratory Tubing

4) Bubble Bottle Assembly



# MOC Question #2

Does 'Bubbling' in the final bubble bottle assembly guarantee CPAP delivery to baby?

True or False

# MOC Answer #2

## **False –**

- 1) Pinched or occluded prongs could cause 'bubbling' without CPAP delivery to baby
- 2) Congested or occluded nasal passages could cause 'bubbling' without CPAP delivery to baby
- 3) Auscultation to listen for CPAP transmission is the standard of care to determine for CPAP delivery

# Building Bubble CPAP



QUESTIONS?